

Development of the First Automated Antibiotic Susceptibility Test System for *Neisseria*

Gonorrhoeae

Sachidevi Puttaswamy*, Mitchell D Newton, Andrew Isakson, Jared McPeake, Rani R Kalam, Ted Hadfield, Roy R Swiger

Acenxion Biosystems, Inc.



ACENXION
BIOSYSTEMS™

Poster #2 Booth #204


BACKGROUND

The Problem

Gonorrhea is the second most commonly reported bacterial infection in the United States. Simultaneously, it represents a critically urgent threat due to its hypermutability leading to pronounced AMR tendencies. **Of the 1.5 million estimated cases annually¹, ~50% are antibiotic resistant.** CDC has repeatedly called for action on NG AMR (in both 2013² and 2019³ AR threat reports), but to date there is limited effort towards development of rapid, automated, phenotypic NG AST: the critical solution to the AMR problem, as emphasized below

2017 Harvard Modeling Study⁴

- Modeled 3 Antibiotics (treatment)
- Prediction: within 15 YRS NONE remain efficacious
- Conclusion: AST on 40% of cases ensures 40YRS efficacy
- * Azithromycin suspended by the CDC 2020



THE SYSTEM

Consumable




- Handheld
- Inexpensive, Room Temp Stable
- Single use, disposable
- 30-well customizable antibiotic cassettes

Instrument



- Benchtop (Small footprint)
- Automated Phenotypic AST results in hrs
- Seamless clinical workflow integration
- Scalable/Stackable (low to high) throughputs

Software




- Advanced Streamlined Software
- Key data acquired and analyzed automatically
- Custom formatted auto report generation
- Easy integration into IT secure communication

RAPID NG AST SYSTEM

- Rapid and actionable results
- Reduce incidences of antibiotic resistance among NG
- Curbs Antibiotic overuse
- Enable directed patient treatment within hours
- Reduce NG prevalence

A STREAMLINED WORKFLOW (Decentralized)

Promotes "Enhanced Surveillance" and "Test-to-Treat"



- Fully automated phenotypic AST results for gonorrhea
- Reduced technician time and labor
- Enables expanded surveillance testing due to ease of use
- Real-time alert isolate information
- Move from empirical to informed patient treatment
- Reduce AMR spread in NG
- Limit sample losses from cryo and shipment exchanges
- Streamlined result, alert isolate and batched reporting across network

CLINICAL WORKFLOW

Collection



Reporting



Analysis



Average results: < 6 hrs

Time and Effort Saved

METHODS

Our technical approach consisted in four (4) key milestones:

- Analysis of current workflows using CDC's GISP protocol
- Demonstration of feasibility using broth dilution AST of *N. gonorrhoeae*.
- Iterative development of consumables and readers controlled by in-house developed API
- Testing of prototype devices demonstrating results that successfully call MIC for QC and AR Isolate Bank challenge strains of NG

Our System:


- Proprietary media allowing for broth culture of NG
- Vital environmental controls for optimum NG growth
- Automated measurement across individually addressed wells
- Customizable reports for MIC and S/I/R interpretation
- Etest® used in parallel for confirmatory testing

Vast prevalence of antimicrobial resistance in NG and its mutagenicity makes NAAT a poor substitute for phenotypic AST

- NAATs are especially adept for fast & sensitive ID needs
- Can be adopted at Point-of-Care
- Inadequate to inform best course of treatment
- Can't yield MICs or give AST profiles
- ESC resistance involves mosaic allele mutations creating major challenges for NAAT

RESULTS

All Graphical Results are Raw, Unprocessed System Data



Results Summary

- Extensive Antibiotic susceptibility testing (AST) of NG strains (10+) including challenge AR isolates from CDC & FDA AR bank across multiple classes of antibiotics
- Yields quantitative results: both MIC and S/I/R
- Average time to call is 6.5 hrs (using simple algorithm)
- Time-to-call Min - 2.33 hrs and Max - 15.67 hrs

Current and Future studies

- Reproducibility testing underway
- Evaluating lyophilized antibiotics
- Azithromycin antibiotic testing currently being investigated
- Test panel expansion to 7 antibiotics across 5-8 concentrations each
- Beginning evaluation of machine learning algorithms to make accurate MIC calls

| Antibiotics | Essential Agreement (EA) vs Agar dilution | Essential Agreement (EA) vs BMD method | Categorical Agreement (CA) |
|---------------|---|--|----------------------------|
| Cefixime | 89% | 100% | 100% |
| Ciprofloxacin | 66.7% | 100% | 100% |
| Tetracycline | 88.9% | 77.8% | 77.8% |



Beta concept Model

SUMMARY

Our system & method is the world's first fully automated, phenotypic AST for *N. gonorrhoeae*. This product addresses the current clinical lab workflows established by the CDC and can be used by public health laboratories at the local, state, and regional levels. Composed of handheld, single-use consumable; small footprint, benchtop instrument; and intuitive and user-friendly graphical-interface; **we show rapid, phenotypic, actionable results with quality-control and challenge strains of NG.**

With the threat of NG developing resistance to all known antibiotics now a reality, **the need to develop targeted, informed treatment regimens has become imperative.** Such regimen is only possible with a testing suite encompassing automated, phenotypic testing which to date has not been shown to be feasible.

We employ patent-pending and proprietary intellectual property and incorporate key elements of existing workflows to demonstrate feasibility for implementation across multiple areas of Public Health interest

For More Information

General Contact:

Roy R Swiger, roy@acxbio.com
Mitchell Newton, mitchell@acxbio.com

Poster Contact:

Sachidevi Puttaswamy, sachi@acxbio.com

Website: <https://acenxionbiosystems.com/>

COME VISIT OUR BOOTH AT #204

References

- <https://www.cdc.gov/std/gonorrhea/drug-resistant/public-health-threat/default.htm>
- <https://www.cdc.gov/drugresistance/pdf/ar-threats-2013-2018.pdf>
- <https://www.cdc.gov/drugresistance/pdf/threats-report-2019-ar-threats-report-2018.pdf>
- Tuite, Ashleigh R., et al. "Impact of rapid susceptibility testing and antibiotic selection strategy on the emergence and spread of antibiotic resistance in gonorrhoea." *The Journal of infectious diseases* 216.9 (2017): 1141-1149.
- Resistance or Elevated Minimum Inhibitory Concentration (MIC) Patterns of *Neisseria gonorrhoeae* Isolates to Antimicrobials, Gonococcal Isolate Surveillance Project (GISP), 2019

Acronyms:

NG - *Neisseria gonorrhoeae*; AMR - Antimicrobial resistance; GISP - Gonococcal Isolate Surveillance Project; MIC - Minimum Inhibitory Concentrations; S/I/R - Susceptible/Intermediate/Resistant; NAAT - Nucleic Acid Amplification Test; ESC - Extended Spectrum Cephalosporins

Acknowledgements / Sources

Funding: Current funding by NIAID: #5R44AI157675

Investors: Centennial Investors, Columbia, Missouri
Other Private Investment

Early Seed funding: CDC contract #75D30118Co2888



CENTENNIAL INVESTORS

THE ACENXION TEAM



Dr. Roy R Swiger
CEO
Dir. Microbiology



Mr. Andrew Isakson
Dir. Engineering



Dr. Sachidevi Puttaswamy
Dir. Systems Integration



Mr. Mitchell Newton
Sr. Scientist



Dr. Ted Hadfield
Sr. Scientist